

Instructions for Use

Version: 1.0.1
Revision date: 5-Aug-25

Helicobacter Pylori Antibody (Anti-HP) Rapid Test Kit

Catalog No.: abx092204

Size: 400 / 2,000 / 10,000 tests

Storage: Store all reagents at 2-30°C. Keep dry.

Application: For qualitative detection of Anti-HP in Human plasma, serum and whole blood.

Introduction and assay principle

Abbexa's Helicobacter Pylori Antibody (Anti-HP) Rapid Test Kit is based on the gold immuno-chromatography assay (GICA) principle. When the concentration of Anti-HP in the sample is more than the detection limit, there is a color change in the detection line and the result is positive. When the concentration of Anti-HP in the sample solution is less than the detection limit, there is no color change in the detection line and the result is negative.

Kit Components

- Test cassettes with pipettes
- Sample Diluent

Material Required But Not Provided

- Timer

Sample preparation

Fresh samples are recommended. Avoid repeated freeze/thaw cycles, bacterial pollution, visible particles; and avoid cloudy, hemolytic, or viscous samples.

- **Serum and Plasma:** Serum and plasma samples should be collected using conventional methods. Test immediately or store between 2 – 8°C for short-term storage (up to 5 days) or -20 °C for long-term storage.
- **Whole blood:** Collect samples in an anticoagulant tube. Test immediately or store at 2-8 °C for up to 24 hours.

Note: Allow samples and Sample Diluent to equilibrate to room temperature before use.

Assay procedure


1. Take a test cassette and lay it flat on a clean table. Using the provided pipette, slowly and vertically add 1 drop (approximately 25 µl) of sample into the sample well on the test cassette. Immediately add 3 drops of sample diluent (approximately 75 µl) of Sample Diluent into the sample well. Avoid foaming.
2. Leave at room temperature for 10 minutes, then analyze the result.

Instructions for Use

Version: 1.0.1
Revision date: 5-Aug-25

Results analysis

- **Positive result:** A colored line is observed in both the control (C) section and the test (T) section.
- **Negative result:** A colored line is observed in the control (C) section but not the test (T) section.
- **Invalid result:** No colored line is observed in the control (C) section.



Notes

1. The test cassettes should be brought to room temperature before use.
2. After opening the aluminum foil, use the test cassette as soon as possible.
3. Samples should be clear with no visible particles, turbidity, or bacterial pollution.
4. Do not mix or re-use the disposable pipettes to avoid cross-contamination.
5. Do not use water, PBS, or similar solutions, as the negative control.
6. Avoid touching the cassette membrane through the sample well or test result window.
7. This kit is for qualitative detection of Anti-HP in Human serum, plasma and whole blood samples. For other sample types, a preliminary experiment is recommended to determine compatibility with this kit.
8. This kit is for research use only and the results are for reference only. A rough estimate of the antibody concentration can be determined by the intensity of the color of the detection line, using the provided antibody titer color card. It is recommended to use this kit in conjunction with another detection method.
9. All waste should be disposed of appropriately. Please note that you may need to follow special waste disposal procedures for infectious samples. Please check local disposal regulations.
10. Do not interpret the results after 15 minutes.

Technical Support

For troubleshooting and technical assistance, please contact us at support@abbexa.com.