

Instructions for Use

Version: 1.0.1

Revision date: 14-Aug-25

6-Monoacetylmorphine (6-MAM) Rapid Test Kit

Catalog No.: abx092237

Size: 100 / 400 / 2000 / 10,000 tests

Storage: Store all reagents at 2°C – 30°C. Keep dry.

Application: For qualitative detection of 6-MAM in human urine.

Introduction and Assay Principle

Abbexa's 6-Monoacetylmorphine (6-MAM) Rapid Test Kit is based on the competitive gold immuno-chromatography assay (GICA) principle. Any 6-MAM present in the samples combines with the colloidal gold particle-labelled 6-MAM antibody, preventing the 6-MAM antibody from binding to immobilized 6-MAM at the detection line. When the concentration of 6-MAM in the sample is more than the detection limit, there is no color change in the detection line and the result is positive. When the concentration of 6-MAM in the sample solution is less than the detection limit, there is a color change in the detection line and the result is negative.

Kit Components

- Test cassettes with pipettes

Material Required But Not Provided

- Specimen collection container
- Timer

A. Sample preparation

- **Urine:** Urine samples should be collected using conventional methods in a clean, dry container and stored between 2 – 8°C for short-term storage (up to 48 hours) or -20 °C for long-term storage. Fresh samples are recommended. Avoid repeated freeze/thaw cycles, bacterial pollution, visible particles; and avoid cloudy, hemolytic, or viscous samples.

B. Assay procedure

- Take a test cassette and lay it flat on a clean table. Using the provided pipette, slowly and vertically add 3 drops (approximately 60 µl) of sample to the sample well on the test cassette. Avoid foaming.
- Leave at room temperature for 5 – 10 minutes, then analyze the result.


Instructions for Use

Version: 1.0.1

Revision date: 14-Aug-25

C. Results analysis

- **Positive result:** : A colored line is observed in the control (C) section but not the test (T) section.
- **Negative result:** A colored line is observed in both the control (C) section and the test (T) section.
- **Invalid result:** No colored line is observed in the control (C) section.



Notes

1. The test cassettes should be brought to room temperature before use.
2. After opening the aluminum foil, use the test cassette as soon as possible.
3. Samples should be clear with no visible particles, turbidity, or bacterial pollution.
4. Do not mix or re-use the disposable pipettes to avoid cross-contamination.
5. Do not use water, PBS, or similar solutions as the negative control.
6. Avoid touching the cassette membrane through the sample well or test result window.
7. This kit is for qualitative detection of 6-MAM in human urine samples. For other sample types, a preliminary experiment is recommended to determine compatibility with this kit. Positive samples can be tested with another method (e.g. HPLC, LC/MS) for quantitative results.
8. This kit is for research use only and the results are for reference only.
9. All waste should be disposed of appropriately. Please note that you may need to follow special waste disposal procedures for infectious samples. Please check local disposal regulations.

Technical Support

For troubleshooting and technical assistance, please contact us at support@abbexa.com.