

**Instructions for Use**

Revision date: 22-Oct-20

**Mammalian Total Protein Extraction Kit****Catalog No.:** abx098853**Size:** 100 ml**Storage:** Store at -20°C for up to one year.

**Introduction:** Abbexa's Mammalian Total Protein Extraction Kit is a fast and efficient way to extract total protein (including cytoplasmic, membrane and nuclear proteins) from mammalian cells and tissues without ultracentrifugation. The extracted proteins can be analyzed by SDS-PAGE, ELISA, WB and other functional assays.

**Kit components**

1. Total Protein Extraction Buffer (TPEB): 100 ml
2. 100X EDTA-free Protease Inhibitor Cocktail: 1 ml

**Materials Required But Not Provided**

1. Phosphate-Buffered Saline (PBS)
2. Phenylmethanesulfonyl Fluoride (PMSF)
3. Distilled water
4. High-precision pipette and sterile pipette tips
5. Centrifuge and centrifuge tubes
6. Timer
7. Ice
8. Vortexer
9. Homogenizer

**Notes**

- Before use, add PMSF (not provided in this kit) to the Total Protein Extraction Buffer (TPEB). Mix thoroughly and aliquot after first use.
- All steps should be carried out on ice or at 2-8 °C.
- A BCA assay is recommended for protein quantification.

**Procedure****A. Cultured cells**

1. Harvest  $0.5-2 \times 10^7$  cells. Centrifuge at  $500 \times g$  for 5 minutes to pellet the cells. Wash the pellet with 1 ml of pre-chilled PBS and centrifuge at  $500 \times g$  for 5 minutes. Carefully discard the supernatant and repeat the wash one more time.
2. Add 1 ml of TPEB to the pellet. Mix thoroughly by vortexing, and allow to stand on ice for 30 minutes, vortexing every 10 minutes.
3. Centrifuge at  $14,000 \times g$  at 4 °C for 10 minutes.
4. Carefully transfer the supernatant (which contains the total protein) to a new 1.5 ml microcentrifuge tube. The isolated total proteins can be analyzed immediately or stored at -80°C.

**Instructions for Use**Revision date: 22-Oct-20

---

**B. Tissues**

1. Cut and mince 20-100 mg of tissue into small pieces. Wash with 1 ml of pre-chilled PBS and centrifuge at 500 × g for 5 minutes. Carefully discard the supernatant and repeat the wash one more time.
2. Add 1 ml of TPEB to the tissue. Mix thoroughly by vortexing. Transfer the suspension to a pre-chilled glass homogenizer and homogenize the tissue with 6-10 strokes.
3. Carefully transfer the suspension to a new 1.5 ml microcentrifuge tube. Allow to stand on ice for 30 minutes, vortexing every 10 minutes.
4. Centrifuge at 14,000 × g at 4 °C for 10 minutes.
5. Carefully transfer the supernatant (which contains the total protein) to a new 1.5 ml microcentrifuge tube. The isolated total proteins can be analyzed immediately or stored at -80°C.

For Reference Only