

Instructions for Use

Revision date: 22-Oct-20

Mammalian Membrane Protein Extraction Kit

Catalog No.: abx098855

Size: 50 rxns

Storage: Store the 100X EDTA-free Protease Inhibitor Cocktail at -20°C and the Membrane Protein Extraction Buffers at 2-8°C for up to one year.

Introduction: Abbexa's Mammalian Total Protein Extraction Kit is a simple and efficient way to extract membrane proteins from mammalian cells and tissues without ultracentrifugation in 70 minutes, with up to 90% efficiency for membrane proteins with 1-2 transmembrane domains. The extracted proteins can be analyzed by SDS-PAGE, ELISA, WB and other functional assays.

Kit components

1. Membrane Protein Extraction Buffer 1 (MPEB1): 50 ml
2. Membrane Protein Extraction Buffer 2 (MPEB2): 7.5 ml
3. Membrane Protein Extraction Buffer 3 (MPEB3): 15 ml
4. 100X EDTA-free Protease Inhibitor Cocktail: 1 ml

Materials Required But Not Provided

1. Phosphate-Buffered Saline (PBS)
2. Phenylmethanesulfonyl Fluoride (PMSF)
3. Distilled water
4. High-precision pipette and sterile pipette tips
5. Centrifuge and centrifuge tubes
6. Timer
7. Ice
8. Vortexer
9. Homogenizer

Notes

- Before use, add Protease Inhibitor Cocktail and PMSF (not provided in this kit) to Membrane Protein Extraction Buffers 1, 2 and 3 (MPEB1, MPEB2 and MPEB3).
- All steps should be carried out on ice or at 2-8 °C.
- A BCA assay is recommended for protein quantification.

Procedure

A. Cultured cells

1. Harvest $0.5-1 \times 10^7$ cells. Wash the cells with 1 ml of pre-chilled PBS and centrifuge at $1000 \times g$ for 3 minutes. Carefully discard the supernatant and repeat the wash one more time.
2. Add 750 μ l of MPEB1 to the cell pellet. Mix thoroughly by vortexing for 15 seconds, and allow to stand on ice for 10 minutes, vortexing every 2 minutes.
3. Centrifuge at $16,000 \times g$ at 2-8 °C for 15 minutes.
4. Carefully transfer the supernatant (which contains cytoplasmic proteins) to a new 1.5 ml microcentrifuge tube. The isolated cytoplasmic proteins can be analyzed immediately or stored at -80°C.
5. Add 150 μ l of MPEB2 to the cell pellet. Resuspend the pellet by vortexing for 15 seconds, and allow to stand on ice for 30

Instructions for Use

Revision date: 22-Oct-20

minutes, vortexing every 5 minutes.

6. Centrifuge at 16,000 × g at 2-8 °C for 15 minutes.
7. Carefully transfer the supernatant (which contains membrane proteins) to a new 1.5 ml microcentrifuge tube. The isolated cytoplasmic proteins can be analyzed immediately or stored at -80°C.

B. Tissues

1. Wash 20-60 mg of tissue with 2 ml of pre-chilled PBS and vortex briefly. Carefully discard the supernatant.
2. Add 1 ml of PBS. Cut and mince the tissue into small pieces. Centrifuge at 5,000 × g for 3 minutes, then carefully discard the supernatant.
3. Add 750 µl of MPEB1 to the tissue. Mix thoroughly by vortexing. Transfer the suspension to a pre-chilled glass homogenizer and homogenize the tissue with 6-10 strokes.
4. Allow to stand on ice for 10 minutes, vortexing every 2 minutes.
5. Centrifuge at 16,000 × g at 2-8 °C for 15 minutes.
6. Carefully transfer the supernatant (which contains cytoplasmic proteins) to a new 1.5 ml microcentrifuge tube. The isolated cytoplasmic proteins can be analyzed immediately or stored at -80°C.
7. Add 150 µl of MPEB2 to the cell pellet. Resuspend the pellet by vortexing for 15 seconds, and allow to stand on ice for 30 minutes, vortexing every 5 minutes.
8. Centrifuge at 16,000 × g at 2-8 °C for 15 minutes.
9. Carefully transfer the supernatant (which contains membrane proteins) to a new 1.5 ml microcentrifuge tube. The isolated cytoplasmic proteins can be analyzed immediately or stored at -80°C.

For Reference Only