Datasheet

Version: 3.0.0 Revision date: 18 Oct 2025



Human Anti-Myeloperoxidase Antibody (Anti-MPO) CLIA Kit

Catalogue No.:abx196448

Human Anti-Myeloperoxidase Antibody (Anti-MPO) Chemiluminescent Immunoassay (CLIA) Kit is a Indirect Chemiluminescent Immunoassay (CLIA) Kit for use with Serum, plasma and other biological fluids.

Target: Anti-Myeloperoxidase Antibody (Anti-MPO)

Research Area: Enzymes and Kinases, Infection Immunity, Cardiovascular Biology, Hepatology

Reactivity: Human

Tested Applications: CLIA

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Storage: Shipped at 4°C. Upon receipt, store the kit according to the storage instruction in the kit's

manual.

Validity: The validity for this kit is at least 6 months. Up to 12 months validity can be provided on

request.

Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5%

within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also

strongly suggested that the whole assay is performed by the same user throughout.

UniProt Primary AC: P05164 (UniProt, ExPASy)

Gene Symbol: MPO

GeneID: <u>4353</u>

OMIM: <u>254600</u>

HGNC: 7218

KEGG: hsa:4353

Ensembl: ENSG00000005381

String: <u>9606.ENSP00000225275</u>

Enzyme Commission Number: EC 1.11.2.2

Datasheet

Version: 3.0.0 Revision date: 18 Oct 2025



Test Range: 0.27 ng/ml - 200 ng/ml

Sensitivity: < 0.14 ng/ml

Standard Form: Lyophilized

Detection Method: Chemiluminescent

Assay Type: Indirect

Assay Data: Quantitative

Sample Type: Serum, plasma and other biological fluids.

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR

THERAPEUTIC PROCEDURES.

The range and sensitivity is subject to change. Please contact us for the latest product information. For accurate results, sample concentrations must be diluted to mid-range of the kit. If you require a specific range, please contact us in advance or write your request in your

order comments.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.