

Fast Coomassie Blue Protein Staining Solution

Catalogue No.: abx098912

Special offer: Free of charge with your order. Limit 1 per order. Offer available while stocks last. Please see the [Rewards and Promotions](#) page for further details.

Fast Coomassie Blue Protein Staining Solution is a ready-to-use Coomassie Brilliant Blue stain for polyacrylamide gels. Staining is achieved within 15 minutes of application, and no fixing, washing or destaining steps are required.

Target: Fast Coomassie Blue Protein Staining Solution

Tested Applications: SDS-PAGE

Form: Liquid

Storage: Store at 4 °C.

Note: This product is for research use only.

Directions for use:

1. Mix the reagent before use.
2. After SDS-PAGE is finished, remove the gel from cassette and submerge it in Fast Protein Staining Solution.
3. Incubate for 15 minutes at room temperature.
4. Remove the dye solution carefully.
5. Image the gel.