

One-Step RT-PCR SuperMix (+dye) (8 kb)

Catalogue No.:abx098028

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. Amplified fragment up to 8kb.

Target: One-Step RT-PCR SuperMix (+dye) (8 kb)

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR

COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

This product is shipped with dry ice.

