Datasheet

Version: 4.0.0 Revision date: 13 Mar 2025



2xRNA Loading Buffer

Catalogue No.:abx098243

This product is used as loading buffer in RNA electrophoresis. It is suitable for denatured or native agarose gel and polyacrylamide gel electrophoresis. This product contains tracking dye- Bromophenol Blue, which can monitor the progress of electrophoresis. To use this product in electrophoresis, EB or other nucleic acid staining dye should be added into RNA sample or gel in advance. Or after gel electrophoresis, stain the gel with EB or other dyes, followed by observation under UV light.

Target: 2xRNA Loading Buffer

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR

COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

