

Beagle dog Plasma (Non-Sterile Sodium EDTA)

Catalogue No.: abx800097

This product is currently in development. The lead time for this product may be several months. Please contact us at info@abbexa.com for an updated lead time before purchasing this product.

Plasma blood fraction from Beagle dog in Sodium EDTA. This product is Non-Sterile.

Target:	Beagle dog Plasma (Non-Sterile Sodium EDTA)
Origin:	Beagle dog
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Form:	Liquid
Purification:	Prepared from beagle dog blood.
Storage:	Store at -20°C. Avoid repeated freeze/thaw cycles.
Buffer:	Sodium EDTA.
Concentration:	> 50 mg/ml
Note:	This product is for research use only. This product is shipped with dry ice.