Datasheet

Version: 3.0.0 Revision date: 25 May 2025



Beagle dog Plasma (Non-Sterile Sodium Heparin)

Catalogue No.:abx800099

This product is currently in development. The lead time for this product may be several months. Please contact us at info@abbexa.com for an updated lead time before purchasing this product.

Plasma blood fraction from Beagle dog in Sodium Heparin. This product is Plasma.

Target: Beagle dog Plasma (Non-Sterile Sodium Heparin)

Origin: Beagle dog

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Form: Liquid

Purification: Prepared from beagle dog blood.

Storage: Store at -20°C. Avoid repeated freeze/thaw cycles.

Buffer: Sodium Heparin.

Concentration: > 50 mg/ml

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC,

THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL

CONSUMPTION.

This product is shipped with dry ice.