

## Beagle dog Plasma (Sterile Sodium EDTA)

Catalogue No.: abx800098

Plasma blood fraction from Beagle dog in Sodium EDTA. This product is Plasma.

**Target:** Beagle dog Plasma (Sterile Sodium EDTA)

**Origin:** Beagle dog

**Recommended dilutions:** Optimal dilutions/concentrations should be determined by the end user.

**Form:** Liquid

**Purification:** Prepared from beagle dog blood. 0.22 µm sterile filtered.

**Storage:** Short term: Store at -20°C. Long term: Store at -80°C. Avoid repeated freeze/thaw cycles.

**Buffer:** Sodium EDTA.

**Note:** THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

This product is shipped with dry ice.

**Directions for use:** Use aseptic techniques to maintain sterility when opening the product. Thaw quickly and mix gently at 37 °C, immediately transfer to ice and use within 24 hrs.