

Human Anti-Glutamate Receptor Ionotropic, NMDA 2A Antibody (Anti-GRIN2A) ELISA Kit

Catalogue No.: abx258783

Human Anti-Glutamate Receptor Ionotropic, NMDA 2A Antibody (Anti-GRIN2A) ELISA Kit is an ELISA Kit for the in vitro quantitative measurement of Human Anti-Glutamate Receptor Ionotropic, NMDA 2A Antibody (Anti-GRIN2A) concentrations in serum, plasma and other biological fluids.

Target:	Anti-Glutamate Receptor Ionotropic, NMDA 2A Antibody (Anti-GRIN2A)
Research Area:	Neuroscience
Reactivity:	Human
Tested Applications:	ELISA
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Storage:	Shipped at 4°C. Upon receipt, store the kit according to the storage instruction in the kit's manual.
Shelf Life:	The validity for this kit is at least 6 months. Up to 12 months validity can be provided on request.
Stability:	The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same user throughout.
UniProt Primary AC:	Q12879 (UniProt , ExPASy)
KEGG:	hsa:2903
String:	9606.ENSP00000379818
Test Range:	3.12 ng/ml - 200 ng/ml
Sensitivity:	< 1.33 ng/ml
Standard Form:	Lyophilized
Detection Method:	Colorimetric
Assay Type:	Indirect

Datasheet

Version: 6.0.0

Revision date: 01 Nov 2025



Assay Data: Quantitative

Sample Type: Serum, plasma and other biological fluids.

Target Type: Antibody

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES.

The range and sensitivity is subject to change. Please contact us for the latest product information.

For accurate results, sample concentrations must be diluted to mid-range of the kit. If you require a specific range, please contact us in advance or write your request in your order comments.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.

For Reference Only