


E-Selectin / CD62E (SELE) Antibody (FITC)

Catalogue No.: abx271276



Fluorescence emission spectra of FITC.

Selectin, Endothelium (SELE) Antibody (FITC) is a Rabbit Polyclonal antibody conjugated to FITC against Selectin, Endothelium (SELE).

Target: E-Selectin / CD62E (SELE)

Clonality: Polyclonal

Reactivity: Pig

Tested Applications: WB, IHC, IF/ICC

Host: Rabbit

Recommended dilutions: WB: 0.5-2 µg/ml, IHC: 5-20 µg/ml, IF/ICC: 5-20 µg/ml. Optimal dilutions/concentrations should be determined by the end user.

Conjugation: FITC

Excitation/Emission: 499/515

Laser Line: 488

Isotype: IgG

Form: Liquid

Purification: Purified by antigen-specific affinity chromatography.

Storage: Aliquot and store at -20°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.

Buffer: 0.01 M PBS, pH 7.4, containing 0.05% Proclin-300, 50% glycerol.

Datasheet

Version: 3.0.0

Revision date: 11 Aug 2025



Note:

THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

For Reference Only