

Coating Buffer

Catalogue No.: abx098970

Coating Buffer formulated to aid capture antibody binding to ELISA plate wells.

Target: Coating Buffer

Tested Applications: ELISA

Form: Lyophilized

Reconstitution: Reconstitute in 220 ml distilled water. If a precipitate is observed, heat in a water bath set at 37 °C until the precipitate dissolves.

Storage: Store lyophilized at room temperature for up to one year. After reconstitution, store at 2-8 °C for up to one week.

Buffer: The exact components are proprietary. Contains Na_2CO_3 and NaHCO_3 .

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

Directions for use: Avoid microbial contamination. Personal protective equipment (PPE), such as lab coats, lab glasses and disposable gloves, are recommended when using this product.