## **Datasheet**

Version: 8.0.0 Revision date: 18 Jun 2025



## **Bimagrumab ELISA Kit**

Catalogue No.:abx395105

Bimagrumab ELISA Kit is a quantitative ELISA kit for detection of Bimagrumab.

Bimagrumab is a monoclonal anti-ActRIIA/IIB antibody that targets ActRIIA and ActRIIB simultaneously. Native ActRIIA and ActRIIB are both mediators of signaling from several TGFβ superfamily members such as activins, myostatin, and GDF11. Hence, pan-blockade by Bimagrumab achieves greater effect than a monoclonal antibody against either receptor alone, i.e. monotherapy with anti-ActRIIA antibody or anti-ActRIIB antibody. Bimagrumab has been explored as a therapy for sporadic inclusion body myositis, primary sarcopenia, and muscle wasting associated with chronic obstructive pulmonary disease.

Target: Bimagrumab

Research Area: Biosimilars

Reactivity: Human

Tested Applications: ELISA

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

**Storage:** Shipped at 4 °C. Upon receipt, store the kit according to the storage instruction in the kit's manual.

Validity: The validity for this kit is 6 months.

**Stability:** The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% within

the expiration date under appropriate storage conditions. To minimize performance fluctuations, operation procedures and lab conditions should be strictly controlled. It is also strongly suggested

that the whole assay is performed by the same user throughout.

**Test Range:** 125 ng/ml - 4000 ng/ml

Sensitivity: < 97.3 ng/ml

Standard Form: Lyophilized

**Detection Method:** Colorimetric

Assay Type: Competitive

Assay Data: Quantitative

Sample Type: Serum and plasma.

## **Datasheet**

Version: 8.0.0 Revision date: 18 Jun 2025



**CAS Number:** 1356922-05-8

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR

THERAPEUTIC PROCEDURES.

The range and sensitivity is subject to change. Please contact us for the latest product information. For accurate results, sample concentrations must be diluted to mid-range of the kit. If you require a specific range, please contact us in advance or write your request in your order comments. Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.