

One-Step RT-PCR SuperMix (+dye) (8 kb)

Catalogue No.: abx098028

One-Step RT-PCR combines the first-strand cDNA synthesis with PCR in the same tube to simplify reaction setup and reduce the possibility of contamination. Only gene-specific primers can be used for One-Step RT-PCR. Amplified fragment up to 8kb.

Target: One-Step RT-PCR SuperMix (+dye) (8 kb)

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

This product is shipped with dry ice.

For Reference Only