

RNase-Free water

Catalogue No.: abx098145

RNase-free Water is prepared with deionized water, which is treated with 0.01% DEPC, followed by high-temperature autoclave. It is suitable for a wide range of RNA-related molecular biology applications.

Target: RNase-Free water

Storage: Store at -20°C. Stable for 12 months from date of receipt.

Note: THIS PRODUCT IS FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC, THERAPEUTIC OR COSMETIC PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION.

For Reference Only